

Product Datasheet

TUBGCP5 Peptide - C-terminal region (orb2001904)

Description	This is a synthetic peptide designed for use in combination with anti-GCP5 Antibody (ARP75842_P050),...
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: NASASSGSDQGPSSRQHTMVSFLKPVLKQIIMAGKSMQLLKNLQCAESTT
MW	112kDa
Uniprot ID	Q96RT8
NCBI	NP_443135
Expiration Date	6 months from date of receipt.