

Product Datasheet

TREX2 Peptide - C-terminal region (orb2001965)

Description	This is a synthetic peptide designed for use in combination with anti-TREX2 Antibody (ARP75407_P050)...
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: ALRGLDRAHSHGTRARGRQGYSLSLGFHRYFRAEPSAAHSAEGDVHTLLL
MW	25kDa
Uniprot ID	Q9BQ50
NCBI	NP_542432
Expiration Date	6 months from date of receipt.