

Product Datasheet

SF3A2 Peptide - N-terminal region (orb2002228)

Description	This is a synthetic peptide designed for use in combination with anti-SF3A2 Antibody (ARP73890_P050)...
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: GGKTGSGGVASSESNRDRRERLRQLALETIDINKDPYFMKNHLGSYECK
MW	51kDa
Uniprot ID	Q15428
NCBI	NP_009096
Expiration Date	6 months from date of receipt.