

## Product Datasheet

### EXOSC5 Peptide - N-terminal region (orb2002382)

<b>Description</b>	This is a synthetic peptide designed for use in combination with anti-EXOS5 Antibody (ARP40939_P050)...
<b>Tested Applications</b>	WB
<b>Preservatives</b>	Lyophilized powder
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Synthetic peptide located within the following region: CSLRHFACEQNLLSRPDGSASFLLQGDTSVLAGVYGPAEVKVSKEIFNKAT
<b>MW</b>	25kDa
<b>Uniprot ID</b>	<a href="#">Q9NQT4</a>
<b>NCBI</b>	<a href="#">NP_064543</a>
<b>Expiration Date</b>	6 months from date of receipt.