
Product Datasheet

eotaxin antibody (orb20025)

Description

Goat polyclonal to eotaxin

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

ELISA

Target

eotaxin

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: No specific signals and low background in Human Colon, Duodenum, Ileum and Lung lysates - our routinely used western blotting protocol does not allow detection of proteins as small as the ligand derived from its precursor or the precursor itself with calculated size of 10.7kDa according to NP_002977.1. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any assay results been reported with other antibodies

Protein Sequence

QDSMKYLDQKSPTPK

Clonality

Polyclonal

MW

10.7

Entrez[6356](#)**Dilution Range**

ELISA: 1:2000

Expiration Date

12 months from date of receipt.