

Product Datasheet

ARHGAP27 Peptide - middle region (orb2002504)

Description	This is a synthetic peptide designed for use in combination with anti-RHG27 Antibody (ARP67368_P050)...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: SQDKQMLYTNHFTQEVPVPAPRSIHKSSQDGDTPAQASPPEEKVPAELD
MW	94kDa
Uniprot ID	Q6ZUM4
Expiration Date	6 months from date of receipt.