

Product Datasheet

COX5B Peptide - middle region (orb2002621)

Catalog Number	orb2002621
Category	Proteins
Description	COX5B Peptide - middle region
Predicted Reactivity	Human
Form/Appearance	Lyophilized powder
Buffer/Preservatives	Lyophilized powder
Protein Sequence	Synthetic peptide located within the following region: LEREIMLAACKGLDPYNVLAPKGASGTREDPNLVPSISNKRIVGCICEED
UniProt ID	P10606
MW	14kDa
Tested applications	WB
Application notes	This is a synthetic peptide designed for use in combination with COX5B Rabbit Polyclonal Antibody (orb331446). It may block above mentioned antibody from binding to its target protein in western blot and/or immunohistochemistry under proper experimental settings.
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
NCBI	NP_001853

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558