

Product Datasheet

FRA10AC1 Peptide - C-terminal region (orb2002872)

Description	This is a synthetic peptide designed for use in combination with anti-FRA10AC1 Antibody (ARP68362_P0...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: HHRRKEIKSKKRKDKTKKDCEESSHKKSRLSSAEEASKKKDKGHSSSKKS
MW	30kDa
Uniprot ID	Q70Z53
Expiration Date	6 months from date of receipt.