

Product Datasheet

ERP29 antibody (orb20031)

Description	Goat polyclonal antibody to ERP29
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Target	ERP29
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx. 28kDa band observed in Human Liver lysates (calculated MW of 29.0kDa according to NP_006808.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	EKNKMSDGGKKEELQ
Clonality	Polyclonal
MW	29.0; 5.48
Entrez	10961
Dilution Range	ELISA: 1:2000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb20031
(0.1ug/ml)
staining of
Human Li...