

## Product Datasheet

### ZBED2 Peptide - middle region (orb2003255)

<b>Description</b>	This is a synthetic peptide designed for use in combination with anti-ZBED2 Antibody (ARP72042_P050)...
<b>Tested Applications</b>	WB
<b>Preservatives</b>	Lyophilized powder
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Synthetic peptide located within the following region: RQVSRGPGVNVGTTALWKHLKSMHREELEKSGHGQAGQRQDPRPHGPQLP
<b>MW</b>	23kDa
<b>Uniprot ID</b>	<a href="#">Q9BTP6</a>
<b>NCBI</b>	<a href="#">NP_078784</a>
<b>Expiration Date</b>	6 months from date of receipt.