

Product Datasheet

PTAR1 Peptide - C-terminal region (orb2003696)

Description	This is a synthetic peptide designed for use in combination with anti-PTAR1 Antibody (ARP71479_P050)...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LISQTVIDSSVMEQNPLRSEPALVPPKDEEAAVSTEEPRINLPHLLEEEV
MW	46kDa
Uniprot ID	Q7Z6K3
NCBI	NP_001093136
Expiration Date	6 months from date of receipt.