

Product Datasheet

OR11H12 Peptide - middle region (orb2003716)

Description	This is a synthetic peptide designed for use in combination with anti-OR5D3 Antibody (ARP71254_P050)...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: YAFIFITVMKMPSTGGRKKAFSTCASHLTAITIFHGTILFLYCVPSKSS
MW	23kDa
Uniprot ID	B2RN74
NCBI	NP_001013372.1
Expiration Date	6 months from date of receipt.