

Product Datasheet

MCCC2 Peptide - C-terminal region (orb2003833)

| | |
|----------------------------|---|
| Description | This is a synthetic peptide designed for use in combination with anti-MCCC2 Antibody (ARP46313_P050)... |
| Tested Applications | WB |
| Preservatives | Lyophilized powder |
| Form/Appearance | Lyophilized powder |
| Storage | Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles. |
| Note | For research use only |
| Protein Sequence | Synthetic peptide located within the following region: RKVVRNLNYQKKLDVTIEPSEEPLFPADELYGIVGANLKR SFDVREVIAR |
| MW | 36kDa |
| Uniprot ID | D6R9R1 |
| Expiration Date | 6 months from date of receipt. |