


Product Datasheet

Der P1 protein (orb20043)

Description	Der P1 protein
Species/Host	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, SDS-PAGE, WB
Target	Der P1
Form/Appearance	50mM Tris, 50mM NaCl, pH 8.5
Concentration	0.5 mg/ml
Storage	Store at -80°C for long or -20°C for 6 months. Aliquot to avoid repeated freezing and thawing.
Note	For research use only
Application notes	Der P1 N-His6 tag protein.
Protein Sequence	<p>RPSSIKTFEE YKKAFNKSYSY TFEDEEAARK NFLESVKYVQ SNGGAINHLS DLSLDEFKNR FLMSAEAFEH LKTQFDLNAE TNACSINGNA PAEIDLRQMR TVTPIRMQGG CGSCWAFSGV AATESAYLAY RNQSLDLAEQ ELVDCASQHG CHGDTIPRGI EYIQHNGVVQ ESYRYVARE QSCRRPNAQR FGISNYCQIY PPNVNKIREA LAQTHSAIAV IIGIKDLDAF RHYDGRTIQ RDNGYQPNYH AVNIVGYSNA QGV DYWIVRN SWDTNWGDNG YGYFAANIDL MMIEEYPYVV IL</p>
Purity	> 90.0% as determined by 12% SDS-PAGE
MW	38.1 kDa
Uniprot ID	P08176
Expiration Date	6 months from date of receipt.



WB analysis of BSA 5ug (Lane1), BSA 10 u...