

Product Datasheet

MRPS5 Peptide - N-terminal region (orb2005010)

Description	This is a synthetic peptide designed for use in combination with anti-MRPS5 Antibody (ARP70888_P050)...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: ETGAGAKKGRGKRTKKKKRKDLNRGQIIGEGRYGFLWPGLNVPLMKNGAV
MW	47kDa
Uniprot ID	P82675
NCBI	NP_114108
Expiration Date	6 months from date of receipt.