

Product Datasheet

DCAF8L2 Peptide - C-terminal region (orb2005378)

Description	This is a synthetic peptide designed for use in combination with anti-DCAF8L2 Antibody (ARP70117_P05...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: VKIWTPTAKAATELTGLKKVIKKNKWERDEDESLHHGSLFDQYMLWFLLRH
MW	71kDa
Uniprot ID	J3KT06
NCBI	NP_001130005
Expiration Date	6 months from date of receipt.