

Product Datasheet

Fahd2a Peptide - C-terminal region (orb2005591)

Description	This is a synthetic peptide designed for use in combination with anti-Fahd2a Antibody (ARP62544_P050...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: VSQFVTLYPGDLLLTGTPPGVGMFRKPPVFLKKGDEVQCEIEELGVIINK
MW	35kDa
Uniprot ID	Q3TC72
NCBI	NP_083905
Expiration Date	6 months from date of receipt.