

Product Datasheet

Nkx3-1 Peptide - middle region (orb2005731)

Description	This is a synthetic peptide designed for use in combination with anti-Nkx3-1 Antibody (ARP37115_P050...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	TPTEPESDAHFFETYLLDCEHNPGDLASAPQVTKQPQKRSRAAFSHTQVIE
MW	27kDa
Uniprot ID	P97436
NCBI	NP_035051
Expiration Date	6 months from date of receipt.