

Product Datasheet

Sfpi1 Peptide - N-terminal region (orb2005896)

Description	This is a synthetic peptide designed for use in combination with anti-Sfpi1 Antibody (ARP37131_P050)...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	MLQACKMEGFSLTAPPSDDLVTYDSELYQRPMHDYYSFVGSDGESHSDDHY
MW	31kDa
Uniprot ID	P17433
NCBI	NP_035485
Expiration Date	6 months from date of receipt.