

Product Datasheet

TAPT1 Peptide - C-terminal region (orb2006085)

Description	This is a synthetic peptide designed for use in combination with anti-TAPT1 Antibody (ARP67839_P050)...
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	KEAKMEEKLSNPPATCTPGKPSSKSNKCKPSQGLSTEENLSASITKQPI
MW	62kDa
Uniprot ID	NM_153365
NCBI	Q6NXT6
Expiration Date	6 months from date of receipt.