

Product Datasheet

PPP3CB Peptide - middle region (orb2006099)

Description	This is a synthetic peptide designed for use in combination with anti-PPP3CB Antibody (ARP57515_P050...
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	MCDLLWSDPSEDFGNEKSQEHFHNTVIRGCSYFYNYPAVCEFLQNNLLS
MW	57kDa
Uniprot ID	NM_017042
NCBI	P20651
Expiration Date	6 months from date of receipt.