



Product Datasheet

MTHFD1 antibody (orb20061)


Description	Goat polyclonal antibody to MTHFD1
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	MTHFD1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:1000.WB: Approx. 110kDa band observed in Human Brain (substantia nigra) lysates (calculated MW of 102kDa according to NP_005947.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	RGDLNDCFIPCTPK
Clonality	Polyclonal
MW	102
Entrez	4522
Dilution Range	ELISA: 1:1000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20061
(3.75µg/ml)
staining of
paraff...



orb20061
(1µg/ml)
staining of
Jurkat (A...



Western blot
analysis of
Human
Substanti...