

Product Datasheet

Clmp Peptide - N-terminal region (orb2006129)

Description	This is a synthetic peptide designed for use in combination with anti-Clmp Antibody (ARP49867_P050),...
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	LWLVYTYVGTLGTHTEIKRVAEEKVTLPCHHQLGLPEKDTLDIEWLLTDN
MW	40kDa
Uniprot ID	NM_173154
NCBI	Q8K1G0
Expiration Date	6 months from date of receipt.