

Product Datasheet

CTTNBP2NL Peptide - N-terminal region (orb2006177)

Description	This is a synthetic peptide designed for use in combination with anti-CTTNBP2NL Antibody (ARP68240_P...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	LEEERQRHAQDTAEGDDVTYMLEKERERLTQQLEFEKSQVKKFEKEQKKL
MW	70kDa
Uniprot ID	Q9P2B4
NCBI	NP_061174
Expiration Date	6 months from date of receipt.