

Product Datasheet

PLEKHG3 Peptide - C-terminal region (orb2006312)

Description	This is a synthetic peptide designed for use in combination with anti-PLEKHG3 Antibody (ARP67751_P05...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	PGGRPSARSPLSPTETFSWPDVRELCSKYASRDEARRAGGGRRGPPVNR
MW	128kDa
NCBI	NP_056364
Expiration Date	6 months from date of receipt.