

Product Datasheet

Wnt11 Peptide - middle region (orb2006533)

Description	This is a synthetic peptide designed for use in combination with anti-Wnt11 Antibody (ARP41258_P050)...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	PGCSCGPVPGPEPPGPGNRWGGCADNLSYGLLMGAKFSDAPMKVKKTSQA
MW	39kDa
Uniprot ID	P48615
NCBI	NP_033545
Expiration Date	6 months from date of receipt.