

Product Datasheet

FEZ2 Peptide - C-terminal region (orb2006715)

Description	This is a synthetic peptide designed for use in combination with anti-FEZ2 Antibody (ARP66800_P050),...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	IEVQNKQKEHKETAKKKKLNKGSSQNGKNERSHMPGTYLTTVIPYEKKN
MW	40kDa
Uniprot ID	Q9UHY8
NCBI	NP_005093
Expiration Date	6 months from date of receipt.