

Product Datasheet

COASY Peptide - C-terminal region (orb2006776)

| | |
|----------------------------|---|
| Description | COASY Peptide - C-terminal region |
| Tested Applications | WB |
| Preservatives | Lyophilized powder |
| Form/Appearance | Lyophilized powder |
| Storage | Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles. |
| Note | For research use only |
| Application notes | This is a synthetic peptide designed for use in combination with COASY Rabbit Polyclonal Antibody (orb586536). It may block above mentioned antibody from binding to its target protein in western blot and/or immunohistochemistry under proper experimental settings. |
| Protein Sequence | ETYRGGMAINRFRENDLEELALYQIQLLKDLRHTENEEDKVSSSSFRQR |
| MW | 65kDa |
| Uniprot ID | Q13057 |
| NCBI | NP_001035997 |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)