

Product Datasheet

PBOV1 Peptide - N-terminal region (orb2006791)

Description	This is a synthetic peptide designed for use in combination with anti-PBOV1 Antibody (ARP66421_P050)...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	RHRLSNFPRLPGILAPETVLLPFCYKVFRKKEKVKRSQKATEFIDYSIEQ
MW	16kDa
Uniprot ID	Q9GZY1
NCBI	NP_067648
Expiration Date	6 months from date of receipt.