

Product Datasheet

PTGFR Peptide - C-terminal region (orb2006828)

Description	This is a synthetic peptide designed for use in combination with anti-PTGFR Antibody (ARP65097_P050)...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	RFYLLLFSLGLLALGVSLLCNAITGITLLRVKFKSQHRQGRSHHLEMV
MW	37kDa
Uniprot ID	P43088
NCBI	NP_000950
Expiration Date	6 months from date of receipt.