

## Product Datasheet

### ATG10 Peptide - C-terminal region (orb2006858)

<b>Description</b>	This is a synthetic peptide designed for use in combination with anti-ATG10 Antibody (ARP59065_P050)...
<b>Tested Applications</b>	WB
<b>Preservatives</b>	Lyophilized powder
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	FFVLHPCKTNEFMTPVLKNSQKINKNVNYITSWLSIVGPPVGLNLPLSYA
<b>MW</b>	25kDa
<b>Uniprot ID</b>	<a href="#">Q9H0Y0</a>
<b>NCBI</b>	<a href="#">NP_001124500</a>
<b>Expiration Date</b>	6 months from date of receipt.