

Product Datasheet

CAPN8 Peptide - middle region (orb2007118)

Description	This is a synthetic peptide designed for use in combination with anti-CAPN8 Antibody (ARP66137_P050)...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	KLIRLRNPWGEVEWSGAWSDDAPEWNHIDPRRKEELDKKVEDGEFWMSLS
MW	77kDa
Uniprot ID	A6NHCO
NCBI	NP_001137434
Expiration Date	6 months from date of receipt.