

Product Datasheet

UBE2V2 Peptide - C-terminal region (orb2007560)

Description	This is a synthetic peptide designed for use in combination with anti-UBE2V2 Antibody (ARP64750_P050...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	GMVDARSIPVLAKWQNSYSIKVVLQELRRLMMSKENMKLPQPPEGQTYNN
MW	16kDa
Uniprot ID	Q15819
NCBI	NP_003341
Expiration Date	6 months from date of receipt.