

Product Datasheet

SRMS Peptide - middle region (orb2008094)

Description	This is a synthetic peptide designed for use in combination with anti-SRMS Antibody (ARP63389_P050),...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	PCMPQKAPRQDVWERPHSEFALGRKLGEGYFGEVWGLWLGSLPVAIKVI
MW	54kDa
Uniprot ID	Q9H3Y6
NCBI	NP_543013
Expiration Date	6 months from date of receipt.