

Product Datasheet

VKORC1 Peptide - C-terminal region (orb2008207)

Description	This is a synthetic peptide designed for use in combination with anti-VKORC1 Antibody (ARP63330_P050...
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	LHVKAARARDRDYRALCDVGTAISCSRVFSSRLPADTLGLCPDAAELPGV
MW	10kDa
Expiration Date	6 months from date of receipt.