

Product Datasheet

GTPBP5 Peptide - C-terminal region (orb2008220)

Description	GTPBP5 Peptide - C-terminal region
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	45kDa
Uniprot ID	Q9H4K7
NCBI	NP_056481
Expiration Date	6 months from date of receipt.