

Product Datasheet

UBE2T Peptide - N-terminal region (orb2008264)

Description	UBE2T,ubiquitin-conjugating enzyme E2T (putative),PIG50; HSPC150; UBE2T,ubiquitin-conjugating enzyme E2 T,ubiquitin-conjugating enzyme E2 T,ubiquitin-protein ligase T,ubiquitin carrier protein T,ubiquitin conjugating enzyme,cell proliferation-inducing gene 50 protein,HSPC150 protein similar to ubiquitin-conjugating enzyme
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	MQRASRLKRELHMLATEPPPGITCWQDKDQMDDLRAQILGGANTPYEKGV
MW	22kDa
Uniprot ID	Q9NPD8
NCBI	NP_054895
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)