

## Product Datasheet

### UBE2T Peptide - N-terminal region (orb2008264)

<b>Description</b>	UBE2T,ubiquitin-conjugating enzyme E2T (putative),PIG50; HSPC150; UBE2T,ubiquitin-conjugating enzyme...
<b>Tested Applications</b>	WB
<b>Preservatives</b>	Lyophilized powder
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	MQRASRLKRELHMLATEPPPGITCWQDKDQMDDLRAQILGGANTPYEKGV
<b>MW</b>	22kDa
<b>Uniprot ID</b>	<a href="#">Q9NPD8</a>
<b>NCBI</b>	<a href="#">NP_054895</a>
<b>Expiration Date</b>	6 months from date of receipt.