

## Product Datasheet

### USP7 Peptide - middle region (orb2008485)

<b>Description</b>	This is a synthetic peptide designed for use in combination with anti-USP7 Antibody, made by Aviva S...
<b>Tested Applications</b>	WB
<b>Preservatives</b>	Lyophilized powder
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	RVLLDNVENKMKGTCVEGTIPKLFGRGKMVSYIQCKEVDYRSDRREDDYDI
<b>MW</b>	128kDa
<b>Uniprot ID</b>	<a href="#">Q93009</a>
<b>NCBI</b>	<a href="#">NP_003461</a>
<b>Expiration Date</b>	6 months from date of receipt.