

Product Datasheet

USP38 Peptide - middle region (orb2008487)

Description	This is a synthetic peptide designed for use in combination with anti-USP38 Antibody, made by Aviva ...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	LVNKDVPQKPGGETTPSVTDLLNYFLAPEILTGDNQYYCENCASLQNAEK
MW	116kDa
Uniprot ID	Q8NB14
NCBI	NP_115946
Expiration Date	6 months from date of receipt.