


Product Datasheet

TARDBP antibody (orb20085)


Description	Goat polyclonal antibody to TARDBP
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	TDP-43
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx. 45kDa band observed in Human lysates of cell line Jurkat (calculated MW of 44.7kDa according to NP_031401.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	HISNAEPKHNSNRQ
Clonality	Polyclonal
MW	44.7
Entrez	23435
Dilution Range	ELISA: 1:16000, WB: 0.3-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb20085 (0.3ug/ml) staining of Jurkat I...



Western blot analysis of Human Cerebellu...



Immunohistochemical staining of Human Pr...