

Product Datasheet

Twf1 Peptide - C-terminal region (orb2008565)

Description	This is a synthetic peptide designed for use in combination with anti-Twf1 Antibody (ARP56488_P050),...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	EKLSKRQLNYVQLEIDIKNETIILANTENTELKDLPKRIPKDSARYHFFL
MW	40kDa
Uniprot ID	Q5RJR2
NCBI	NP_001008521
Expiration Date	6 months from date of receipt.