

Product Datasheet

Tnik Peptide - middle region (orb2008647)

Description	This is a synthetic peptide designed for use in combination with anti-Tnik Antibody (ARP48771_P050),...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	DSGAEYGIGSSTKASFTPFVDPRVYQTSPTDEDEEDESSAAALFTSELL
MW	145kDa
Uniprot ID	P83510
NCBI	NP_001156480
Expiration Date	6 months from date of receipt.