

Product Datasheet

SEC23IP Peptide - N-terminal region (orb2009073)

Description	This is a synthetic peptide designed for use in combination with anti-SEC23IP Antibody(ARP52295_P050...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	VPFIPVTQASASPASLLLPGEDSTDVGEEDSFLGQTSIHTSAPQTFSYFS
MW	111kDa
Uniprot ID	Q9Y6Y8
NCBI	NP_009121
Expiration Date	6 months from date of receipt.