

Product Datasheet

RGL3 Peptide - middle region (orb2009222)

Description	This is a synthetic peptide designed for use in combination with anti-RGL3 Antibody (ARP56251_P050),...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	LDWKVLFCHVHLVSIFLSALNSSANPIIYFFVGSFRQRQNRQNLKLVLQR
MW	36kDa
Uniprot ID	Q3MIN7
NCBI	NP_001030300
Expiration Date	6 months from date of receipt.