

Product Datasheet

Rapgef2 Peptide - middle region (orb2009278)

Description	This is a synthetic peptide designed for use in combination with anti-Rapgef2 Antibody (ARP42770_P05...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	SILPQKPYNDIGIGSQDDSIIVGLRQTKHIPAALPVSGTLSSSNPDLLQS
MW	150kDa
Uniprot ID	E9QNQ4
NCBI	NP_001093094
Expiration Date	6 months from date of receipt.