

Product Datasheet

PRDM16 Peptide - middle region (orb2009420)

| | |
|----------------------------|---|
| Description | This is a synthetic peptide designed for use in combination with anti-PRDM16 Antibody(ARP33260_P050)... |
| Tested Applications | WB |
| Preservatives | Lyophilized powder |
| Form/Appearance | Lyophilized powder |
| Storage | Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles. |
| Note | For research use only |
| Protein Sequence | LNHTQDAKLPSPLGNPALPLVSAVSNSSQGTAAAGPEEKFESRLEDSCV |
| MW | 140kDa |
| Uniprot ID | Q9HAZ2 |
| NCBI | NP_071397 |
| Expiration Date | 6 months from date of receipt. |