

Product Datasheet

Ppp1r1a Peptide - N-terminal region (orb2009445)

Description	Ppp1r1a Peptide - N-terminal region
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	19kDa
Uniprot ID	Q9ERT9
NCBI	NP_067366
Expiration Date	6 months from date of receipt.