

Product Datasheet

PKN2 Peptide - C-terminal region (orb2009511)

Description	This is a synthetic peptide designed for use in combination with anti-PKN2 Antibody (ARP61894_P050),...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	TVFDIQNDRNSILPKSQSEYKPDTPQSGLEYSIGIQELEDRRSQRFQFNL
MW	108kDa
Uniprot ID	Q16513
NCBI	NP_006247
Expiration Date	6 months from date of receipt.