

Product Datasheet

Oxsr1 Peptide - middle region (orb2009635)

Description	This is a synthetic peptide designed for use in combination with anti-Oxsr1 Antibody (ARP48661_P050)...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	RAPTISERSKKVRRVPGSSGRLHKTEDGGWEWSDDDEFDEESEEGRAAISQ
MW	52kDa
Uniprot ID	Q6P9R2
NCBI	NP_598746
Expiration Date	6 months from date of receipt.